

A^{Pro}
ADVANTAGE[®]

by **NDC**

Sterile

Cauteries

Pro Advantage high temperature cauteries are used for procedures such as pinpoint hemostasis, evacuation of a subungual hematoma, controlling diffuse bleeders, vasectomies, and sculpting woven grafts. Pro Advantage low temperature cauteries are used primarily for ophthalmology.

- Battery operated
- Lightweight, disposable
- Provide consistent temperature
- Locking safety cap prevents accidental operation
- Four year shelf life
- Individually packaged, sterile
- High temperature available in three tip options
- Low temperature available with Fine tip



Item #	Description
P212000	Ophthalmic Fine Tip Cautery, Low Temperature
P212050	Loop Tip Cautery, High Temperature
P212100	Fine Tip Cautery, High Temperature
P212150	Vasectomy Tip Cautery, High Temperature

Distributed By: