

Fecal Occult Blood Test Devices



Rapid, convenient, and qualitative methods for asymptomatic occult bleeding. Available in guaiac slide – with both color enhanced sensitivity and readability.

- On board positive and negative monitors
- Includes enhanced color sensitivity and readability
- CLIA waived and FDA 510(k) cleared to market



CPT Code: 82270

Item #	Description
P080018	Fecal Occult Blood Lab Pack – Guaiac Slides <i>Includes: 100 Single Slides Tests, 2x10mL Developers, Applicators. 100 per box</i>
P080028	Fecal Occult Developer Solution – for Guaiac Slides <i>Includes: 8x10mL Developers. 8 per box</i>

Distributed By: