

Whole Blood / Serum / Plasma *H. pylori* Test Device

**Pro
ADVANTAGE**

by **NDC**

A one step, qualitative membrane strip based immunoassay for the detection of *H. pylori* IgG antibodies in whole blood, serum or plasma.

H. pylori is one of the most common bacterial infections in humans causing active and chronic gastritis, non-ulcer dyspepsia and gastric/duodenal ulcer formation. The Pro Advantage *H. pylori* Test Device – Whole Blood/Serum/Plasma offers a simple and quick test format which assists the physician in making a treatment decision during the same office visit.

- 88% sensitivity; 89% specificity
- Accurate results in 10 minutes
- Reduces the risk of developing more serious illnesses
- Lessens the likelihood of complications
- CLIA waived for whole blood and FDA 510(k) cleared to market
- Fast and easy – no biopsy or histology expertise required



**For Whole Blood,
Serum/Plasma:
Moderately Complex**
CPT Code:
86318QW Whole Blood,
86677 Serum/Plasma

Item #	Description
P080095	H. pylori Test Device – Whole Blood/Serum/Plasma <i>Includes: 20 individually packaged devices, 20 disposable sample droppers, 20 disposable heparinized capillary tubes, 1 dispensing bulb, 1 sample buffer, 1 procedure card, 1 package insert. 20 tests per box</i>

Distributed By: