

## Bulb & Piston

# Irrigation Syringes


### Bulb Irrigation Syringe

Contoured finger grips and pliable bulb with textured surface provide stable grip and allow maximum control of drawn or expelled fluids. Clear and easy-to-read graduations are in 5 cc (mL) increments, up to 60 cc. Comes with tip protector. Sterile, packaged in poly pouch. Latex-Free.



### Piston Irrigation Syringes

Ideal for irrigating feeding tubes and adding liquid medications to an existing gastrostomy tubes. Available with thumb ring or flat top. Finger rests help provide control. Precision crafted rubber O-ring ensures a positive seal every time the syringe is drawn. 60cc capacity barrel with graduations in oz and cc. Comes with small tip adapter and tip protector. Packaged in a resealable IV pole bag. Latex-Free.



Item #	Description	Box
	Bulb Irrigation Syringe	
<b>P250600</b>	<b>60cc Bulb Irrigation Syringe, Catheter Tip with tip protector</b>	<b>50</b>
	Piston Irrigation Syringe	
<b>P250610</b>	<b>60cc Flat Top Piston Irrigation Syringe, Catheter Tip, with small tip adapter</b>	<b>30</b>
<b>P250615</b>	<b>60cc Thumb Control Ring Piston Irrigation Syringe, Catheter Tip, with small tip adapter</b>	<b>30</b>

**Distributed By:**