

Plain and Frosted

Microscope Slides and Cover Glass

Microscope Slides

- Guaranteed pre-cleaned
- Premium quality glass resistant to corrosion and fogging
- Tissue between each slide prevents sticking
- Precision polished edges on all 4 sides
- Rounded corners prevent injury
- Measures 25mm x 75mm x 1mm and available in plain or frosted



Cover Glass

- Highest quality white glass
- Free from ragged edges, oily film, striae and fingerprints
- Desiccant in each 1-oz. box prevents fogging and sticking
- Measures 22mm x 22mm and available in #1 and #2 thickness



Item #	Description	Quantity
P460025	Plain Microscope Slides	72/box
P460125	Frosted Microscope Slides	72/box
P460122	Cover Glass (0.13 to 0.17mm thick)	1 ounce/box
P460222	Cover Glass (0.17 to 0.25mm thick)	1 ounce/box