

Whole Blood / Serum / Plasma

Infectious Mononucleosis Test Device

Pro
ADVANTAGE[®]

by **NDC**

Rapid chromatographic immunoassay technology identifies heterophile antibodies in whole blood, serum or plasma to aid in the diagnosis of Infectious Mononucleosis.

- Qualitative membrane strip based immunoassay for detection of IM heterophile antibodies
- All inclusive 20 test kit
- Simple to perform using whole blood (from venipuncture or fingerstick), serum or plasma
- Reliable results in 5 minutes
- CLIA waived for whole blood and FDA 510(k) cleared to market



**For Whole Blood
Serum/Plasma:
Moderately Complex**
CPT Code:
86308QW Whole Blood,
86308 Serum/Plasma

Item #	Description
P080016	Infectious Mononucleosis Test Device <i>Includes: 20 individually packaged devices, 20 disposable sample droppers, 20 disposable heparinized capillary tubes, 1 dispensing bulb, 1 sample buffer, 1 procedure card, 1 package insert. 20 tests per box</i>

Distributed By: