

Push-Button

Aluminum Crutches

Features push-button adjustment and patient height scale printed on the crutches for fast, accurate fitting. Includes underarms pads, handgrips and tips. Wing nut adjustment for hand grips.

- Latex-free
- Push-button allows for quick and easy adjustment
- Constructed of high-quality anodized aluminum
- Fully assembled with pad, handgrip and tip in place
- Lightweight, strong and easy-to-clean
- Weight capacity: 350 lbs.



Item #	Description	Crutch Height
P501031	Pediatric	31.5" – 37.5" (fits patient 4' – 4'6")
P501037	Youth	37" – 45" (fits patient 4'6" – 5'2")
P501045	Adult	45" – 53" (fits patient 5'2" – 5'10")
P501053	Tall Adult	53" – 61" (fits patient 5'10" – 6'6")