



**PRO
ADVANTAGE**



Orthopedic Rehabilitation Products Guide

QUALITY PRODUCTS FOR HEALTHCARE PROFESSIONALS

*Manufactured for NDC, Inc.
407 New Sanford Road / La Vergne, TN 37086
www.ProAdvantagebyNDC.com*

ORTHOPEDIC SOFT GOODS



Cervical Collars

- Features anatomically contoured foam with soft cotton stockinette covering for maximum patient comfort
- Indicated for cervical facet syndrome, discogenic pain, degenerative disc disease and muscle sprains or strains of cervical region
- Manufactured from medium firm density foam; 2"W contact closure
- Sizing: Fits 16"-24"

P662050 2 1/2", Universal

P663050 3 1/2", Universal



Arm Slings

- One piece design for maximum comfort, slide buckle closure for ease of application; thumb loop and neck pad provide added comfort
- Indicated for injuries to the arm, shoulder, wrist or hand and as a cast support to maintain arm in neutral position
- Manufactured from a durable cotton material; 1"W shoulder strap with pad

P664010 Small, 15" x 6", Cotton

P664020 Medium, 17" x 7", Cotton

P664030 Large, 19" x 8", Cotton

P664040 X-Large, 21" x 9", Cotton



Shoulder Immobilizers

- Features shoulder and waist straps both with slide buckle closure; thumb loop and neck pad provide added comfort
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery and other shoulder related injuries
- Manufactured from durable cotton material; 1"W shoulder and waist straps

P665010 Small, 15" x 6", Cotton

P665020 Medium, 17" x 7", Cotton

P665030 Large, 19" x 8", Cotton

P665040 X-Large, 21" x 9", Cotton



Clavicle Straps

- Features heavy-duty felt center point pad with hook and loop closure for optimal comfort
- Indicated for clavicle fractures, Grade II/III acromioclavicular separations, sternoclavicular dislocations, postural control
- Constructed of breathable, open cell stockinette covered foam with contact closure

P667010 X-Small, 17" - 24"

P667020 Small, 25" - 31"

P667030 Medium, 27" - 33"

P667040 Large, 29" - 35"

P667050 X-Large, 31" - 37"

ORTHOPEDIC SOFT GOODS



Sling and Swathe

- Features cradle-style sling and swathe with contact closure
- Indicated for immobilization following rotator cuff repairs, S/P glenohumeral acromioplasty surgery or other shoulder-related injuries
- Manufactured from seamless, durable, cloth-laminated foam material

P666000 Universal, Up to 53" Chest Circumference



Lacer Wrist Orthosis

- Anatomically designed and made from high quality materials for enhanced durability and patient compliance; universal sizing reduces inventory
- Indicated for Carpal Tunnel Syndrome, tendonitis, sprains or strains, soft tissue injury, wrist contracture and post-cast removal
- Constructed of breathable, cottony-soft, moisture-wicking liner and high performance brushed nylon laminate exterior; malleable dorsal and palmar aluminum stays
- Sizing: Fits 5"-9"; measurement taken from wrist circumference

P668010 Universal, Right

P668015 Universal, Left

Lacer Wrist with Abducted Thumb Orthosis

- Adjustable cinch closure features two separate pull straps for a secure fit with minimal effort; malleable thumb stay provides extension and abduction control; contoured palmar stay supports the wrist and is easily modified
- Indicated for deQuervain's Syndrome, tendonitis, ulnar collateral ligament injury, CMC arthritis, scaphoid injuries, strains or sprains and post-cast removal
- Constructed of breathable, cottony-soft, moisture-wicking liner and high performance brushed nylon laminate exterior; four malleable aluminum stays
- Sizing: Fits 5"-9"; measurement taken from wrist circumference

P668020 Universal, Right w/ Abducted Thumb

P668025 Universal, Left w/ Abducted Thumb



Compression Knee Immobilizer

- Features contoured posterior stays, movable medial/lateral stays and elastic compression straps at all closures for optimal immobilization
- Indicated for post-op immobilization and mild to severe knee injuries where full extension of the leg is desired
- Constructed of breathable foam laminate; concave aluminum stays; premium quality elastic straps with contact closure

P669016 16" Universal, 28" Thigh Circumference

P669018 18" Universal, 28" Thigh Circumference

P669020 20" Universal, 29" Thigh Circumference

P669022 22" Universal, 29" Thigh Circumference

*Quality Products for
Healthcare Professionals*

*Manufactured for NDC, Inc.
407 New Sanford Road
La Vergne, TN 37086*

© 2016 NDC, Inc.
www.ProAdvantagebyNDC.com