

Disposable

Specimen Containers

Disposable sample collection containers have standard 4 ounce fill capacity and graduations in oz and cc (ml). Leak resistant. Available in sterile and non-sterile.

- Latex Free
- Non-Toxic



Item #	Description	Box
Sterile P250400	4 oz. Specimen Container with screw on lid & label, packaged; Individually in Poly Pouch	100
P250410	4 oz. Specimen Container with screw on lid & Tamper Evident Label	100
Non-Sterile P250415	4 oz. Specimen Container with screw on lid	500

Distributed By: